

BBa_B0015

Part Summary

Double terminator consisting of BBa_B0010 and BBa_B0012. This is the most commonly used terminator in iGEM. It seems to be reliable. Note, however, that Part:BBa_B0014 is a better part for forward and reverse termination.

Sequence

```
ccaggcatcaaataaaacgaaaggctcagtcgaaagactgggcctttcgttttatctgttgttcggtgaacgctc  
tctactagagtcacactggctcaccttcgggtgggcctttctgcgtttata
```

Part Type

Basic Part

Related Parts


B0010, B0012, B0014

Pigeon Image



B0010 B0012

Plasmid Map



Designer Information

Author(s)

Reshma Shetty

Date

July 17, 2003

Affiliation

Stanford University

Contact

endy@stanford.edu

Design Details

Type	Double terminator
Vector	pSB1C3
Design Components	B0010, B0012

Assembly Information

Assembly Method(s)	BioBricks
Chassis	Escherichia coli
Assembly RFC	10
Scars	yes