

BBa_B0015


Part Summary	Double terminator consisting of BBa_B0010 and BBa_B0012. This is the most commonly used terminator in iGEM. It seems to be reliable. Note, however, that Part:BBa_B0014 is a better part for forward and reverse termination.
Sequence	ccaggcatcaaataaaacgaaaggctcagtcgaaagactgggcctttcgttttatctgttgttgcggtgaacgctc tctactagagtcacactggctcaccttcgggtgggcctttctgcgtttata
Part Type	Basic Part
Related Parts	B0010, B0012, B0014

Pigeon Image



B0010 B0012

Plasmid Map



Designer Information

Author(s)	Reshma Shetty
Date	July 17, 2003
Affiliation	Stanford University
Contact	endy@stanford.edu

Design Details

Type	Double terminator
Vector	pSB1C3
Design Components	B0010, B0012

Assembly Information

Assembly Method(s)	BioBricks
Chassis	Escherichia coli
Assembly RFC	10
Scars	yes